



Webster-pak[®] Sealing Mouse[™] Makes sealing a breeze

The **Sealing Mouse** is the ultimate solution for sealing your Webster-pak foil to the blister. With its ergonomic design, it delivers fast, accurate, and precise sealing, ensuring a perfect air and moisture-resistant seal every time. Users who rely on the Sealing Mouse daily absolutely love it – sealing has never been easier!

KEY FEATURES

- Three heat settings for efficient and reliable foil sealing.
- A **PINK ergonomic handle**, designed for comfortable use, applies heat from the warming plate to the back of the foil.
- In just seconds, gentle pressure and heat securely bond the foil to the blister.
- The **BLUE Pressure Handle** is then used to instantly cool the foil against the blister, ensuring a flawless seal every time.

Experience effortless sealing with the Sealing Mouse – precision and ease, perfected.







Webster-pak Sealing Mouse Makes sealing a breeze

Compatible with:







Webster-pak 28 Low Temp



Unit Dose 7[®]



Unit Dose 7 Webstercard™

↘ Why you'll love the NEW Sealing Mouse

- The quickest and easiest way to seal foil onto a Webster-pak blister.
- Updated ergonomic handle designed for improved safety and usability, suitable for both left and right-handed users.
- Rapid heating to the chosen sealing temperature.
- New temperature dial with three preset options for greater consistency.
- Robust and durable design with high-usage, long-lasting components.
- User-tested to ensure reliability and effectiveness.
- Backlit on/off switch illuminates when the device is active, enhancing awareness of the heated plate.

To learn more, contact us at 1800 244 358

webstercare.com.au 1800 244 358 | info@webstercare.com.au



© (01/2012) Webstercare (ABN 29 974 510 173) – 2025. Webstercare, W Device, Webster-pak & Unit Dose 7 are Registered Trademarks of Manrex Pty Ltd t/as Webstercare, whilst Sealing Mouse & Webstercard are Trademarks of Manrex Pty Ltd t/as Webstercare. V1

Webster-pak Community